



View Test Reports... Three easy steps...

1



Open the CNS Vital Signs Application... CLICK MENU ...CLICK VIEW REPORT

2



SELECT the Patient Record and CLICK the 'View Report' Button

3

The Report will present in your Browser or in a PDF form... CLICK Print and print the report for use during the clinical visit or import the report into your EMR - electronic medical record.

The screenshot shows a detailed test report for a patient. The report includes patient information such as name, age, gender, and date of birth. It also contains a table of test results with columns for 'Patient Problem', 'Assessment', 'Findings', and 'Recommendations'. The table is color-coded with green, yellow, and red cells to indicate the status of various tests.

How do I View Reports?